



## New SAA Test Products

Plasvacc is excited to announce our partnership with VMRD as an official distributor of feline diagnostics.

### **vmrd** DIAGNOSTIC TESTS

#### **Feline Serum Amyloid A Test**

The VMRD Feline SAA test measures serum amyloid A (SAA) in feline serum or plasma. SAA is an inflammatory biomarker that helps to diagnose infection and monitor treatment efficacy in multiple species, including cats. It is more sensitive and reliable than body temperature and WBC count for detecting acute infections, as it is highly dynamic in response to systemic inflammation. As infection resolves, SAA decreases rapidly.

This 10-minute lateral flow test can easily be run in-clinic, with numerical results provided by the handheld VMRD Reader. Coupled with other routine tests, measuring feline SAA in-clinic gives veterinarians an invaluable aid for differential diagnosis, early detection of infection, evaluating prognosis, and monitoring clinical condition.

#### **VMRD Reader**

Reader for use with VMRD point-of-care tests that provides a objective result in whole blood and serum samples. Compatible with both equine and feline Serum Amyloid A tests.

#### **HOW TO ORDER**

All products are available for purchase direct from Plasvacc.  
Contact us at:

**P** +61 7 5463 7600

**E** [AUS.Office@Plasvacc.com](mailto:AUS.Office@Plasvacc.com)

