



## New SAA Test Products

Plasvacc is excited to announce our partnership with VMRD as an official distributor of equine diagnostics.

### **vmrd** DIAGNOSTIC TESTS

#### **Equine Serum Amyloid A Test**

Early identification of infection can be challenging, as the clinical signs are often subtle. The VMRD Equine SAA test is an indispensable tool to aid veterinarians in this challenge, as serum amyloid A (SAA) is rapidly responsive to clinical changes. It is virtually undetectable in normal animals but increases with acute inflammation more rapidly, dramatically, and predictably than fibrinogen or WBC count. SAA starts dropping as soon as disease begins to resolve, allowing objective monitoring of clinical condition and treatment efficacy.

VMRD SAA for Equine is a 10-minute horse-side lateral flow test designed for use in field or hospital settings. Test pouches are set up for use with whole blood, with numerical results provided by the VMRD Reader.

#### **VMRD Reader**

Reader for use with VMRD point-of-care tests that provides a quantitative result in whole blood and serum samples. Compatible with both equine and feline Serum Amyloid A tests.

#### **HOW TO ORDER**

All products are available for purchase direct from Plasvacc. Contact us at:

**P** +61 7 5463 7600

**E** [plasmail@plasvacc.com](mailto:plasmail@plasvacc.com)

