



# Product Information

Alvedia's innovative immunochromatographic technology is your source for easy and reliable canine and feline blood typing.

Based in Lyon, France, Alvedia provides a full range of canine and feline diagnostics that are easy to handle, easy to read, and provide reliable results.

## CANINE DIAGNOSTIC KITS

### Canine Blood Typing

These easy-to-use in-clinic tests are designed to determine the major immunoreactive blood types of dogs. The tests utilise immunochromatography and monoclonal antibodies to DEA1.1 (dog). These tests can be performed with as little as one mL of blood, with results obtained within three minutes.

Available individually (QuickTest) or in a pack of 20 (LabTest).



## FELINE DIAGNOSTIC KITS

### Feline Blood Typing

These easy-to-use in-clinic tests are designed to determine the major immunoreactive blood types of cats. The tests utilise immunochromatography and the A and B blood group antigens (cat). These tests can be performed with as little as one mL of blood with results obtained within three minutes.

Available individually (QuickTest) or in a pack of 20 (LabTest).

## HOW TO ORDER

Alvedia diagnostic test kits can be purchased directly from Plasvacc by email or phone.

Visit our website to learn more about Alvedia's Canine & Feline Diagnostic Kits:



SCAN HERE FOR CANINE BLOOD TYPING PRODUCT DETAILS



SCAN HERE FOR FELINE BLOOD TYPING PRODUCT DETAILS

